

With the SNaP[®] Wound Care System, perfection can be found in the details.

A closer look reveals the innovative features that create powerful wound therapy, without any power.

Silent Cartridge

- Clinically proven negative pressure offered at 125 mmHg, 100 mmHg and 75 mmHg
- Small, silent, lightweight design disappears under clothes
- Proprietary spring mechanism generates consistent, even levels of pressure
- No noisy motors, electronics or batteries to interfere with daily living
- Visual indicator displays when cartridge is full or discharged
- BioLock™ technology to gel exudate for improved exudate containment

Advanced Wound Dressing

- Disposable components with off-the-shelf convenience
- Proprietary hydrocolloid dressing and anti-microbial gauze offer periwound protection and easy removal
- Fully integrated micro-port enables flexibility and a tight bending radius for wounds located in difficult areas
- Cut-to-length tubing with one-way flow valve for improved safety

Comfort Strap

- Strap enables device to be conveniently and discreetly worn for maximum discretion

