


- Every patient should have an invisible friend.



# Say hello to SNaP.<sup>®</sup>

My partner in healing.




The SNaP<sup>®</sup> Wound Care System combines the portability of advanced wound dressings with the proven efficacy of negative pressure therapy in a discreet design that won't get noticed.

Active healing that's out of sight.<sup>™</sup>



# People won't even notice its elegant design.

The SNaP® Wound Care System stays hidden from the outside world while providing patients optimal negative pressure wound therapy throughout their day.



- Small, silent, lightweight design disappears under clothes
- No noisy motors, electronics or batteries to interfere with daily living
- Visual indicator when cartridge is full or discharged
- Disposable components with off-the-shelf convenience

	SNaP®	KCI ActiV.A.C.®*	KCI V.A.C. Via™†
Length	15.1 cm	15.2 cm	17.2 cm
Width	5.5 cm	19.3 cm	9.9 cm
Depth	2.7 cm	6.4 cm	6.1 cm
Weight	66 g	1080 g	320 g

\*Information obtained at [www.kci1.com/KCI1/activatherapyunit](http://www.kci1.com/KCI1/activatherapyunit)

†Information obtained at [www.vacvia.com](http://www.vacvia.com)



# You'll find its engineering simply brilliant.

With the SNaP® Wound Care System, perfection can be found in the details. A closer look reveals the innovative features that create powerful wound therapy, without any power.

- Proprietary spring mechanism generates consistent, even levels of pressure
- Cut-to-length tubing with one-way flow valve for improved safety
- Proprietary hydrocolloid dressing and anti-microbial gauze<sup>†</sup> offer periwound protection and easy removal

† Kendall Kerlix™ AMD™

## Success revealed.

The SNaP® Wound Care System may be discreet, but when it comes to performance, you can see the results for yourself.



- Clinically proven negative pressure offered at 125 mmHg, 100 mmHg and 75 mmHg
- Demonstrated non-inferiority in wound healing outcomes for patients completing at least 4 weeks of therapy<sup>1</sup>
- Wounds not typically treated with negative pressure therapy may heal faster from earlier intervention enabled by the SNaP® System's unique design

## Imagine therapy without boundaries.

The SNaP® Wound Care System makes a visible difference in patients' quality of life, thereby increasing compliance and improving outcomes.


- The SNaP® System interferes significantly less with overall activity, sleep and social interactions than the V.A.C.® Therapy System<sup>1</sup>
- Unrestricted mobility enables patients to move freely while actively treating their wounds with non-invasive negative pressure
- Continuous therapy over a 24 hour period may improve outcomes

### Exit Survey Data<sup>1</sup>

PERCENTAGE OF PATIENTS WHO...	SNaP®	KCI ActiV.A.C.®
<i>Felt that their NPWT device was "never" or "rarely" noticed in social situations</i>	100%	25%
<i>Felt that noise from their NPWT device "never" bothered them</i>	100%	33%
<i>Felt that their NPWT device did not interfere with their daily activities</i>	92%	42%
<i>"Agreed" or "strongly agreed" that their NPWT device was comfortable to wear</i>	85%	33%


# Savings will be visible.

The SNaP® Wound Care System earns its name by saving you considerable time over both modern wound dressings and intrusive negative pressure devices.



- The SNaP® System could save substantial treatment costs in addition to allowing patients greater freedom and mobility
- Requires only two dressing changes per week, freeing 33% of nurses time<sup>1,2</sup>
- Quicker to apply than a powered negative pressure device [8.58 minutes instead of 20.98]<sup>1</sup>
- Off-the-shelf availability eliminates the rental procurement process

**Cost Savings with the SNaP® Wound Care System vs. Other Wound Therapies<sup>2</sup>**



# Do you have a patient who needs an invisible friend?

Introduce more patients to active wound healing with the SNaP® Wound Care System, the quiet companion that keeps attention on results.

- The SNaP® Wound Care System is indicated for patients who would benefit from a suction device particularly as the device may promote wound healing. The SNaP® Wound Care System is indicated for the removal of small amounts of exudate from the following types of wounds:
  - Chronic
  - Acute
  - Traumatic
  - Subacute
  - Dehisced
  - Ulcers
  - Surgically Closed Incisions
- Available in 3 different pressure settings to meet a variety of clinical needs

SNaP® Wound Care Cartridge		
Catalog Number	Pressure	Quantity
SNP-A-125/5	125 mmHg	5
SNP-A-100/5	100 mmHg	5
SNP-A-075/5	75 mmHg	5

SNaP® Wound Care Dressing Kit		
Catalog Number	Size	Quantity
SKT-A-15x15/10	15 cm x 15 cm	10

SNaP® Wound Care Strap		
Catalog Number	Size	Quantity
STP-A-S/5	Small	5
STP-A-M/5	Medium	5
STP-A-L/5	Large	5

# Say hello to SNaP.<sup>®</sup>

My partner in healing.

Help your patients get their life back on track with the SNaP<sup>®</sup> Wound Care System. It's the best kept secret in wound therapy for people who want active healing.




Active healing that's out of sight.<sup>™</sup>

#### REFERENCES

1 Armstrong, D. G., W. A. Marston, et al. "Comparison of Negative Pressure Wound Therapy with the SNaP<sup>®</sup> Wound Care System vs. V.A.C.<sup>®</sup> Therapy System for the Treatment of Chronic Lower Extremity Ulcers: A Multicenter Randomized Controlled Trial." *Wound Rep Reg* 2011; 19: 173-180.

2 Hutton, D. W. and Sheehan, P. (2011), Comparative effectiveness of the SNaP<sup>®</sup> Wound Care System. *International Wound Journal*, 8: 196-205. doi: 10.1111/j.1742-481X.2011.00775.x



SNAP™

Active healing that's out of sight.™

[www.spiracur.com](http://www.spiracur.com)

---

Spiracur Inc.  
1180 Bordeaux Drive  
Sunnyvale, CA 94089  
Toll Free 1.877.774.7228

Spiracur, Spiracur logo and SNaP are registered trademarks of Spiracur Inc. SNaP & Design and *Active healing that's out of sight* are trademarks of Spiracur Inc. The SNaP® Wound Care System is protected by one or more U.S. patents, with other U.S. and certain foreign patents pending. ©2011 Spiracur Inc. All rights reserved.

L21798 rev. 110419

